

# **Measurement of magneto-optical properties of materials with a Kerr spectrometer**

*M.S. Creative Component  
Iowa State University  
Department of Physics & Astro.  
Kurt Jensen July 1995*

# Dielectric & conductivity tensors

For an otherwise optically isotropic magneto-optic material:

$$\vec{\tilde{e}} = \begin{pmatrix} \tilde{\mathbf{e}} & -i\tilde{Q}\tilde{\mathbf{e}} & 0 \\ i\tilde{Q}\tilde{\mathbf{e}} & \tilde{\mathbf{e}} & 0 \\ 0 & 0 & \tilde{\mathbf{e}} \end{pmatrix} \quad \vec{\tilde{s}} = \begin{pmatrix} \tilde{\mathbf{s}}_{xx} & -\tilde{\mathbf{s}}_{xy} & 0 \\ \tilde{\mathbf{s}}_{xy} & \tilde{\mathbf{s}}_{xx} & 0 \\ 0 & 0 & \tilde{\mathbf{s}}_{xx} \end{pmatrix}$$

$Q$  is the magneto-optic Voigt parameter.

These tensors are related:

$$\vec{\tilde{e}} = \mathbf{e}_0 + i \frac{4p\vec{\tilde{s}}}{W}$$

# Complex Kerr angle

$$\begin{aligned}\tilde{\mathbf{f}}_K &\equiv \mathbf{q}_K + i\mathbf{e}_K = \frac{\tilde{E}_y}{\tilde{E}_x} = -\frac{i\tilde{n}\tilde{Q}}{1-\tilde{n}^2} = -\frac{\tilde{\mathbf{e}}_{xy}}{(1-\tilde{\mathbf{e}})\sqrt{\tilde{\mathbf{e}}}} \\ &= -\frac{4pi}{w} \frac{\tilde{\mathbf{s}}_{xy}}{\tilde{n}(1-\tilde{n}^2)}\end{aligned}$$

In a circular basis, this can be expressed as:

$$\tilde{\mathbf{f}}_K = i \frac{\tilde{n}_+ - \tilde{n}_-}{\tilde{n}_+ \tilde{n}_- - 1} = i \frac{\tilde{r}_+ - \tilde{r}_-}{\tilde{r}_+ + \tilde{r}_-}$$

↑  
Via Fresnel Eqn.

# Complex reflectivity

The complex Fresnel reflectivity at normal incidence is:

$$\tilde{r}_\pm = \tilde{r}_x \pm i\tilde{r}_y = \frac{1 - \tilde{n}_\pm}{1 + \tilde{n}_\pm}$$

where the subscripts denote helicity in the usual sense.  
The Kerr parameters in terms of the reflectivity are:

$$q_K = \frac{\Delta_- - \Delta_+}{2} \quad ; \quad e_K = \frac{r_- - r_+}{r_- + r_+}$$

where

$$\tilde{r}_\pm \equiv r_\pm e^{i\Delta_\pm}$$

# The magnitude of the Kerr effect is proportional to magnetic field

$$n_- - n_+ = n(w + w_L) - n(w - w_L) = 2 \frac{dn}{dw} w_L$$

$$w_L = \frac{eH}{2mc}$$

is the Larmor frequency

so  $f_K = i \frac{\tilde{n}_+ - \tilde{n}_-}{\tilde{n}_+ \tilde{n}_- - 1} \propto H$

# Instrument model

The Jones vector for the detected light is modeled as follows:

$$T = ASMP$$

$$P = \begin{pmatrix} \cos y \\ \sin y \end{pmatrix}$$

$P$  is the Jones vector for the light transmitted by the polarizer.  $y$  is the polarizer angle.

$$M = \begin{pmatrix} 1 & 0 \\ 0 & e^{id} \end{pmatrix}$$

$M$  is the Jones matrix for the modulator.  
 $d$  is the modulator retardation.

$$S = \begin{pmatrix} 1 & 1 \\ -i & i \end{pmatrix} \begin{pmatrix} r_+ e^{i\Delta_+} & 0 \\ 0 & r_- e^{i\Delta_-} \end{pmatrix} \begin{pmatrix} 1 & i \\ 1 & -i \end{pmatrix}$$

$S$  is the sample reflectivity.

$$A = \begin{pmatrix} 1 & 0 \\ 0 & 0 \end{pmatrix} \begin{pmatrix} \cos j & \sin j \\ -\sin j & \cos j \end{pmatrix}$$

$A$  is the Jones matrix for the analyzer.  $\phi$  is the analyzer angle.

# Instrument model

The transmitted intensity is the square magnitude of  $T$ :

$$I = T^* T$$

We modulate the birefringence in M as  $\mathbf{d} = \mathbf{d}_0 \sin \omega t$ . Multiplying the matrices out yields the following for components of  $I$  at DC,  $w$ , and  $2w$ .

$$I_{DC} = 2R[1 + \cos(2y) \cos(2j + 2q_K) + J_0(\mathbf{d}_0) \sin(2j + 2q_K)]$$

$$I_w = -2\Delta R J_1(\mathbf{d}_0) \quad I_{2w} = 4R J_2(\mathbf{d}_0)$$

Where  $2R = r_+^2 + r_-^2 ; \Delta R = r_+^2 - r_-^2$

$$r_+ r_- \approx R ; q_K = \frac{\Delta_- - \Delta_+}{2}$$

# Error propagation

The Kerr signal is proportional to the AC component normalized to the DC level:

$$f = C \frac{V_{AC}}{V_{DC}} \quad ; \quad V_{DC} = V_L - V_D$$

We measure the variance of each measured voltage to find the total variance in the Kerr signal.

$$(df)^2 = C^2 \left( \frac{\cancel{f}}{\cancel{V}_{AC}} \right)^2 (dV_{AC})^2 + \left( \frac{\cancel{f}}{\cancel{V}_L} \right)^2 (dV_L)^2 + \left( \frac{\cancel{f}}{\cancel{V}_D} \right)^2 (dV_D)^2$$

This gives

$$(df)^2 = C^2 \frac{(dV_{AC})^2}{V_{DC}^2} + \frac{V_{AC}^2}{V_{DC}^4} [(dV_L)^2 + (dV_D)^2]$$